
EAR & FOREHEAD THERMOMETER 101

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1 INTRODUCTION

The Ear & Forehead thermometer uses infrared technology to measure infrared energy emitted from the forehead or eardrum and surrounding tissue, and then converts it into temperature value.

2 ABOUT BODY TEMPERATURE

Normal human body temperature fluctuates throughout the day; generally, lowest in the morning, highest in the evening.

Different parts of human body also have different temperatures. Therefore, the standards of measuring oral, axillary, rectal and ear temperatures are different. The human body temperature is controlled and adjusted by the hypothalamus in the brain. As the veins flow through the hypothalamus, they also flow through the ear-drum in the mid ear. Thus, measuring the ear temperature (the temperature of the ear drum membrane) can best reflect the body temperature.

Normal body temperature is a range, not a fixed value.

When using an IR ear thermometer to take temperature, the normal body temperature range is 35,8°C - 38,0°C. Also, people's normal temperature ranges tend to decrease with their age.

However, the normal range is also different from person to person. To determine your and your family members' normal temperature ranges is very important. It can be done easily by using our Infrared Thermometer. Please practice to take temperatures on yourself and your healthy family members to determine their normal temperature ranges.

3 SAFETY NOTES

WARNING

- Please read carefully through the following information concerning safety and proper use. Make yourself familiar with all the functions of the equipment. Be careful to keep these advice notes and if necessary pass them on to a third party.
- This thermometer is intended for household use only. The use of this thermometer is not intended as a substitute for physician's consultation.
- Don't touch or blow the infrared sensor. A polluted infrared sensor may cause inaccuracy.
- Clean a polluted infrared sensor with a soft wet cloth in a gentle manner. Cleaning with toilet paper or paper towel may scratch the infrared sensor, causing inaccuracy.
- Keep the unit out of children's reach. A self-measurement forced by children may cause ear damage. For accidental swallowing of a battery, please contact a doctor at once.
- Measurement must be suspended until the ear is warmed up after outdoor activities to avoid lower results.

- After a shower or exercise, the body temperature will rise. Please take a rest before taking temperature.
- For proper measurement, condition the thermometer for 30 minutes in measurement site till room temperature.
- Earwax or obstructions will affect the accuracy of the thermometer.

NOTE

- Stop using the product in occurrence of any pain. It is likely to damage the ear hole.
- It is not recommended to use the product on users suffering ear disease including otitis external and lymphantitis. Affected part is likely to deteriorate.
- Don't use the product in a wet ear hole after swimming or bathing. It is possible to damage the ear hole.

PROHIBITED

- Don't throw batteries into fire. Batteries may explode in fire

4 CARE AND MAINTENANCE

- Do not shock, vibrate or drop the thermometer.
- Do not allow any liquid in the thermometer.
- Do not disassemble, repair or remodel the thermometer.
- Keep out of reach of children.
- Clean the main body of the thermometer with a soft and dry cloth.
- Do not use benzene, gasoline, thinner or any other cleaner containing abrasive agents.
- The tip of the probe is the most delicate part of the thermometer, accurate readings from the thermometer can be obtained only if the tip of the probe is clean and intact.
- To clean the tip of the probe, hold the thermometer downwards so that no liquid can enter. Use a cotton swab moistened with alcohol to wipe gently and carefully on the surface of the tip of the probe. Allow an hour drying time before operating the thermometer again.


5 BUTTONS

See folded cover page

- | | |
|----------------------------------|--------------------|
| 1. Ear-/Forehead mode detector | 7. LED flash light |
| 2. Reset switch | 8. LIGHT - button |
| 3. Battery compartment (2 x AAA) | 9. LCD display |
| 4. °C/°F switch | 10. START - button |
| 5. OK - button | 11. Forehead probe |
| 6. SET - button | |

6 BATTERY INSTALLATION

1. Remove the battery cover.
2. Install 2 AAA (1,5V) non-rechargeable batteries following the polarity as shown.
3. Reposition the battery cover.

Replace the battery when '  ' is displayed. Take out the battery if you do not use it for a long time.

You have to dispose of the batteries in an environment friendly manner according to your country regulations.

7 HOW TO USE

- It is recommended that you take three temperatures under the following circumstances.
 - Infants in the first 3 months of life.
 - Children under three years old who have a condition such as compromised immune system and for whom the presence or absence of fever is critical.
 - When you are first learning to use the thermometer until you are comfortable with the technique and are gaining consistent readings
- Wait about 30 seconds before measuring the same person to avoid excessive cooling of the skin.
- If the results are different, use the highest reading.
- Clean the probe cover to ensure accurate measurement.
- As with any type of thermometer, slight temperature variations (+/- 0,2 - 0,3°C) can happen, if consecutive measurements are taken.

7.1 EAR TEMPERATURE

- Pull off the Forehead probe from the unit.
- Press the **START** - button (10) to turn the power on. The display will show all symbols.

Before inserting the probe:

For children < 1 year: Pull the ear down and back.

For children > 1 year to adult: Pull the ear up and back.

- While holding the ear, insert the probe to seal the ear canal.
- Press the **START** - button (10) once. You will hear a short beep.
- When you hear a double beep, the measurement is finished.
- Take out the thermometer from the ear canal.
- The result is displayed.

7.2 FOREHEAD TEMPERATURE

- Put on the Forehead probe.
- Gently position the probe onto the center of the forehead.
- Press the **START** - button (10) once. You will hear a short beep.
- When you hear a double beep, the measurement is finished.
- The result is displayed.

At completion of measurement, the unit automatically returns to time mode after 30 seconds. The display will rotate 180°C to make it readable when positioned in the cradle.

The display will automatically turn off completely, 1 minute after the last key pressed (if sleep mode is enabled).

7.3 MEMORY

The result of each measurement will be stored automatically into the memory. If the memory is full (20 records), the oldest measurement result will be removed.

To recall the memory:

- Press the **START** - button (10) to turn the power on.
- Press the **SET** - button (6) successively to scroll through the memory. Each result is displayed with date and time.

7.4 DISPLAY ERROR MESSAGES

'HI'	The measured temperature is higher than 42,9°C.
'LO'	The measured temperature lower than 32°C
'ErrE'	The ambient temperature is not within the operation specifications.
'ErrU'	Does not measure after pressing 'Measure'. Follow the steps mentioned in the user guide.
'ErrH'	An error is detected during self-check. Remove the batteries for 10 seconds and try again.
'ErrS'	Sensor error. Do not use the thermometer under extreme environment.
'ErrR'	Hardware problem. Return the unit to the service center.
'ErrP'	Hardware problem. Return the unit to the service center.

8 ADVANCED MODE

8.1 TIME SETTINGS

- Turn the unit with the probe facing down.
- Power on the unit by pressing the **START** - button (10) once. You will hear a short beep. Time mode is displayed.
- Press and hold the **SET** - button (6). The 24hr/12hr system is blinking.
- Press the **SET** - button (6) to switch between 24hr and 12hr system.
- Press the **OK** - button (5) to confirm. The hour digits are blinking.
- Press the **SET** - button (6) repeatedly to increase the digit.
- Press the **OK** - button (5) to confirm. The minutes digits are blinking.
- Repeat previous steps for setting the minutes, year, month and day.
- After setting the day, the unit will enter sleep mode setting.

8.2 SLEEP MODE - ROOM TEMPERATURE

When sleep mode is not activated, the display is continuously ON.
The display will toggle between time and room temperature.

After setting the hour, sleep mode setting is displayed:

- Press the **SET** - button (6) to switch between **ON** or **OFF**.
- Press the **OK** - button (5) to confirm.

8.3 BACKLIGHT AND FLASHLIGHT

- Press the **LIGHT** - button (8) to turn the display backlight on. The backlight will automatically turn off after 8 seconds.
- Press and hold the **LIGHT** - button (8) to activate the flashlight. The flashlight will automatically turn off after 5 seconds.

9 TEMPERATURE ENTITY

- Power on the unit by pressing the **START** - button (10) once. You will hear a short beep. Time mode is displayed.
- Open the battery compartment.
- Use a small tool to push the **°C/°F Switch** (4) on the back side of the unit.
- A single beep confirms that you have set the temperature to °F.
- A double beep confirms that you have set the temperature to °C.

10 RESET

- Power on the unit by pressing the **START** - button (10) once. You will hear a short beep. Time mode is displayed.
- Open the battery compartment.
- Use a small tool to push the **RESET Switch** (7) on the back side of the unit.
- All segment will be displayed and all stored information and settings are reset.

11 DISPOSAL OF THE DEVICE (ENVIRONMENT)

At the end of the product lifecycle, you should not throw this product into the normal household garbage but bring the product to a collection point for the recycling of electrical and electronic equipments. The symbol on the product, user guide and/or box indicate this.

Some of the product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities in case you need more information on the collection points in your area.

12 FEATURES

Model	HV-T46
Temperature Sensor type	Infrared
Measuring temperature range	32°C ~ 42,9°C
Room temperature	5°C ~ 59,9°C
Display accuracy	0,1°C
Measurement accuracy	35,5°C ~ 42,0°C (+/- 0,2°C) 32,0°C ~ 35,4°C (+/- 0,3°C) 42,1°C ~ 42,9°C (+/- 0,3°C)
	Room temperature +/- 1°C
Operating ambient temperature range	15°C ~ 40°C
Operating ambient humidity range	30 ~ 80% RH
Storage environment temperature range	-10°C ~ +55°C
Storage environment humidity range	30 ~ 85% RH
Memory	Last 20 measurements
Battery Type	2 x AAA Non- rechargeable - DC 3V
Automatic switch-off	1 Minute after last measurement
Power consumption	< 30 mA
Weight	110g (with battery)
Classification 	• Type B • Class IIa • Internally powered equipment • IPX0 • Not suitable for use in presence of flammable anaesthetic mixture with air or with oxygen or nitrous oxide • Continuous operation with short-time loading

13 CE MARK

The CE symbol on the unit, user guide and giftbox indicates that the unit complies with the essential requirements of the medical directive 93/42/EEC.

The declaration of conformity can be found on:

<http://www.topcom.net/support/cedclarations.php>

14 TOPCOM WARRANTY

14.1 WARRANTY PERIOD

The Topcom units have a 24-month warranty period. The warranty period starts on the day the new unit is purchased. Consumables or defects causing a negligible effect on operation or value of the equipment are not covered.

The warranty has to be proven by presentation of the original purchase receipt, on which the date of purchase and the unit-model are indicated.

14.2 WARRANTY HANDLING

A faulty unit needs to be returned to a Topcom service centre including a valid purchase note.

If the unit develops a fault during the warranty period, Topcom or its officially appointed service centre will repair any defects caused by material or manufacturing faults free of charge.

Topcom will at its discretion fulfil its warranty obligations by either repairing or exchanging the faulty units or parts of the faulty units. In case of replacement, colour and model can be different from the original purchased unit.

The initial purchase date shall determine the start of the warranty period. The warranty period is not extended if the unit is exchanged or repaired by Topcom or its appointed service centres.

14.3 WARRANTY EXCLUSIONS

Damage or defects caused by incorrect treatment or operation and damage resulting from use of non-original parts or accessories not recommended by Topcom are not covered by the warranty.

The warranty does not cover damage caused by outside factors, such as lightning, water and fire, nor any damage caused during transportation.

No warranty can be claimed if the serial number on the units has been changed, removed or rendered illegible.

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